

....mct injector

mesocarboxytherapy
by mesoestetic®



mesoestetic®

Technology

mesotherapy

"Little, rarely and in the right place".

Dr. Michel Pistor

Mesotherapy is a medical treatment that involves the intradermal administration of active substances through multiple micro-injections for correcting aesthetic concerns (skin aging, wrinkles, stretch marks, cellulite, alopecia, among others) or as an additional treatment for painful conditions (rheumatic disease, sports medicine, etc.).

Its effectiveness is based on three main concepts:

- **physical stimulation** produced by the needle (mechanical stimulus);
- **pharmacological effect** produced by the solution or a specific blend of active substances;
- **micro-doses** of the solution which activate a greater number of the receptors in the skin.

It is a non-obtrusive procedure requiring no preparation time with the advantage that patients can immediately resume their daily activities.

carboxytherapy


Carboxytherapy is a non-surgical procedure that uses a fine needle to inject medical grade CO₂ into the subcutaneous tissue.

The CO₂ is quickly and widely distributed throughout bordering tissue, employing a powerful vasodilatory effect at the microcirculatory level, increasing blood flow in muscle and adipose tissue and releasing serotonin, bradykinin, histamine and catecholamine. These substances activate the beta-adrenoceptors, stimulating lipolysis of bordering tissue.

In the presence of CO₂, the body will tend to alleviate the instability in local oxygenation. Oxyhemoglobin increases localised oxygen supply and this produces a vasodilatory effect and enhances tissue oxygenation (*Bohr effect*).

The toxins released during the process are eliminated through the lymphatic system, resulting in improved skin tone.

Mesotherapy is the application technique with the highest penetration level of active ingredients:



mesotherapy + carboxytherapy = mesocarboxytherapy

Combined treatments with mct injector

Treatment areas:



wrinkles on the forehead
and glabella

periorbital area
periocular area
preauricular area



neck

hands
neckline



facial contours

nasolabial fold

lip area



arms

back

abdomen
lumbar region
flanks

buttocks

thighs

knees

In a single device, the **mct injector** brings together two complementary techniques: **mesotherapy** and **carboxytherapy**. Thanks to the combination of these two techniques in the same session or in alternate session, the physical therapist can use a complete tool to perform multiple facial and body treatments:

Face, neck, neckline, scalp, back of hands:

- Normalises skin pigmentation;
- Provides additional moisturising;
- Smoothes wrinkles and expression lines;
- Plumps the tissues;
- Redefines the facial contour;
- Improves blood microflow, especially in the scalp;
- Stimulates local microflow, especially in the eye contour.

Arms, back, thorax, abdomen, periumbilical area, flanks, thighs:

- Effectively treats oedematous-fibrosclerotic panniculopathy (cellulite);
- Stimulates the lymphatic system in order to eliminate excess liquids and toxins;
- Activates the lipolysis process;
- Restores firmness to drooping tissue.

Several mesotherapy techniques

mct injector incorporates 3 programs for the administration of active substances through different injection techniques used in traditional mesotherapy.

Meso

For **burst**, **intradermal nappage** or **point by point** injections. **mct injector** automatically adjusts the depth and speed of the needle in order to perform steady and precise movements.

Expert

Program specially designed for **point by point** injections or to perform **intradermal papules**. With this program, professionals can freely adjust the depth (from 1 to 13 mm) and the dose (from 0.01 ml and up).

Filler

Injection program enabling the placement of **implants** or **fillers**. It also enables **tracer mesotherapy** to be carried out, regulating the volume injected, the speed of retrograde injection and the injection depth.



Meso

Expert

Filler

Latest-generation of carboxytherapy

mct injector incorporates a portable carbon dioxide injection system in order to perform **carboxytherapy** treatments comfortably and easily. **Carboxytherapy** plays a coadjuvant role of **mesotherapy** in the treatment of cellulite and adipose tissue but it can also complement practically all treatments requiring improved microcirculation and increased tissue oxygenation (under-eye puffiness, double chins, slow scarring, vascular disease, etc.).

Total control

Handling control over the depth, angle and volume of each injection.

2 injection methods

Manual for the face and periocular area and automatic for the body areas.

Intuitive

Automatic adjustment of gas flow based on the size of the needle.

Painless

Progressive injection with simultaneous micromassage.



Carboxytherapy



Maximum comfort for the doctor

Handling

Ambidextrous design that is both lightweight and **ergonomic**. All accessories are portable (including CO₂ cartridges).

Simplicity

Treatments can be programmed in a matter of seconds thanks to the **touch wheel** and high definition **colour screen**.

Autonomy

Operates equally with a battery (Li-ion) or connected to a power source (90 V - 240 V).

Flexibility


Compatible with all kinds of needles and syringes.
Does not require the use of an exclusive kit.

Visibility

Directional lighting of treated area.

Reliability

Healthcare product with CE marking.



Maximum comfort for the patient

mct injector allows totally painless treatments for the patient thanks to:

- high-speed needle projection;
- no vibratory effects;
- simultaneous micromassage or following injections to improve the diffusion of active substances.

....mct injector

..... **Mesocarboxytherapy:** combination of two techniques with proven effectiveness.

..... **Intuitive, versatile** and **precise** technology.

..... Multiple **body** and **facial treatments**.

..... **Painless** and **safe** treatments.

..... **Medical device** with CE marking.



mesoestetic®

Mesoestetic Pharma Group, s.l.
C/ Tecnologia, 25
08840 - Viladecans (Barcelona)
Tel. +34 902 26 20 31
www.mesoestetic.com
made in Spain



Medical devices are controlled by independent official agencies which certify their quality through the CE mark. Require this safety and efficacy guarantee.

mesoestetic® guarantee: The best reference.
Over 25 years providing solutions and innovating in the medical cosmetics sector.
Spanish company located in Barcelona present in over 60 countries.



Scan the QR code
to see the main
mesotherapy and
carboxytherapy
techniques.