




**solutions for
women's intimate
treatment**

▽ intimate medical solutions

Gyneco-aesthetics is a medical sub-speciality dealing with the treatment of cosmetic, regenerative and functional concerns of the vulvar-vaginal area.

in mesoestetic® we work with specialists in gynaecology to improve the health and intimate wellness of women. Our research has resulted in solutions to restore functionality, comfort and appearance of the intimate area, ultimately improving the quality of life.

Oestrogens: the key in intimate aging

Oestrogen reduction is responsible for specific histological changes in womens genitourinary area (trigone, urethra, vagina, vulva, vestibule and labia) which affect the optimum functionality and elasticity of tissues.



Age is a critical factor in oestrogen activity levels with symptoms highest during climacteric, however there are other external factors that can cause losses in the quality of tissues and functionality of the area.⁽¹⁾⁽²⁾

Risk factors:

Dependent on oestrogen levels

> CLIMACTERIC

- > anti-oestrogen medicines
- > ovary removal



Independent of oestrogen levels



Treatments
(surgery, radio/chemotherapy)



Tight clothes



Pregnancy/Labour



Practising certain sports
(e.g., horse riding, biking)




Weight changes




Laser/wax epilation

Histopathological changes in women's intimate area

Vaginal area



Vulvar area



▽ intimate medical solutions

Solutions with hyaluronic acid

Products with specific rheological profiles to moisturise and restore women's intimate area adapting to the viscoelastic properties of each application area.



1 syringe X 1 ml


mesofiller[®] sensitive

INTIMATE FILLER MEDICAL SOLUTIONS

Cross-linked hyaluronic acid gel for the treatment of vaginal atrophy, dryness, dyspareunia and sexual dysfunction.

Concentration [cross-linked HA]	20 mg/ml
G' Elastic module	196 Pa
G'' Viscous module	25 Pa
Application areas	vaginal vestibule, introitus, vaginal walls, clitoris, G spot
How does it work?	It improves functionality of the vaginal mucosa, promotes regeneration of atrophic tissue and increases the sensitivity of erogenous zones.

Vaginal mucosa



Improves hydration of the vaginal mucosa. Enhances transport and cell trophic function. Stimulates angiogenesis. Activates collagen and elastin synthesis. Increases epithelial thickness. Restores the volume of erogenous zones enhancing friction during SR (clitoris and G spot).⁽³⁾



2 syringes X 1 ml

mesofiller® sensitive plus

INTIMATE FILLER MEDICAL SOLUTIONS

Cross-linked hyaluronic acid gel for the treatment and correction of vulvar ptosis and hypertrophy in the external genital area.

Concentration [cross-linked HA]	25 mg/ml
G' Elastic module	305 Pa
G'' Viscous module	40 Pa
Application areas	labia majora
How does it work?	It fills and restores tissue volume in the labia majora to preserve their protective function and improve their appearance.



single-dose 12 x 5 ml

hydragel intimate

MOISTURISING SOLUTIONS

Moisturising gel for women's intimate area.

Moisturiser

Hyaluronic acid at different molecular weights to return hydration and juiciness to the tissue.

Rejuvenating

It contains encapsulated retinol that activates cell differentiation.

Biorebalancer


It helps maintain the natural balance of the area due to the presence of a biopolysaccharide.

▽ intimate medical solutions

mesofiller® sensitive: vaginal health parameters⁽⁴⁾ [HA 20 ml]


Gloria Bachmann Vaginal Health Index (VHI):

The results show a significant improvement in the functional parameters of the vagina (**elasticity, flow volume, epithelial integrity, pH, hydration**)




Vaginal maturation index (VMI):

The result shows the **re-establishment of the number of** superficial and intermediate **cells**, which means an improvement of vulvovaginal atrophy. Measured by cytology.





Sexual function (FSFI):

A reduction in the FSFI index shows a significant improvement in the parameters of **sexual arousal, orgasm, satisfaction**.




Nugent index:

A reduction of this index shows a re-establishment of **normal values of non-pathological bacteria**. Measured by cytology.



Visual analogue scale (VAS):

Measurement of pain reduction during the gynaecological examination after the treatment.



Evaluation of satisfaction degree


Patient satisfaction questionnaire

mesofiller® sensitive plus: clinical results⁽⁵⁾

[HA 25 ml]



CASE 1:

Age: 44



CASE 2:

Age: 40




CASE 3:

Age: 48




CASE 4:

Age: 34



▽ intimate medical solutions

Depigmentation solutions



[meso] intimate system*
Six active ingredients with complementary, synergistic action encapsulated in liposomes for maximum efficacy.


dermamelan® intimate

PIGMENT CONTROL SOLUTIONS

Depigmenting medical method developed specifically to act on the conditions responsible for skin hyperpigmentation in women's intimate area.

Application areas

- > genital-perineal area
- > mons pubis
- > inner thighs and groins



protocol and products included in the pack



4 ml
dermamelan® intimate peeling



50 ml
post-peel neutralizing spray



7 ml
dermamelan® intimate mask



7 ml
dermamelan® intimate occlusive film



TREATMENT AT THE CLINIC

intensive depigmentation

STEP 1

After cleaning the area and protecting sensitive parts* with **mesolips protector**, apply **dermamelan® intimate peeling** and leave on for 7 minutes. After this time, neutralise and spray with **post-peel neutralizing spray**.



exposure time 7 min

STEP 2

Apply a large layer of 3-4 mm in thickness of **dermamelan® intimate mask** over the area to be treated, then apply **dermamelan® intimate occlusive film** over the treated area. Once the film is applied, the patient may return home. Leave on for 2 hours.




2 h

*mucosa, labia minora and their external area.

Mechanism of action

dermamelan® intimate exerts a corrective action on the **accumulated melanin** and a regulating action on **inflammation**, a characteristic component of hyperpigmentation, as well as **anti-aging action**, visibly improving skin quality.



10 ml
post-procedure crystal
fiber intimate mask



50 ml
dermamelan® intimate home
depigmenting gel cream



2 h



HOME TREATMENT

Regulation of hyperpigmentation

Apply **post-procedure crystal fiber intimate mask** over the area treated and leave on for 10-15 minutes. After 48 hours apply evenly one or two pumps of **dermamelan® intimate home depigmenting gel cream** all over the area. Use **dermamelan® intimate home depigmenting gel cream** every night for the time recommended by your medical professional.



Removal of dermamelan® intimate occlusive film and dermamelan® intimate mask


Remove the mask with plenty of lukewarm water and using an intimate cleansing gel with an adequate pH.

▽ intimate medical solutions


dermamelan intimate®: studies *in vitro*

Evaluation of the efficacy of [meso] intimate system®

Inhibition of **melanin production** in a melanocyte culture.




Inhibition of the **production of IL-6**, the main inflammatory marker, evidencing the regulatory capacity of the inflammatory component associated with hyperpigmentation.




Evaluation of the efficacy of home depigmenting gel cream

> DEPIGMENTING ACTION



Reduction of tyrosinase enzyme⁽¹²⁾




Inhibition of melanosomes transfer⁽¹³⁾



Reduction of melanin formed⁽¹⁴⁾

> BARRIER FUNCTION




Transepidermal water loss (TEWL*)⁽¹⁵⁾

studies *in vitro*

study *in vivo*

dermamelan® intimate: clinical results⁽¹⁶⁾

CASE 1



CASE 2



CASE 3



CASE 4



▽ intimate medical solutions

Solutions for women's intimate treatment that provide **tangible results and a high degree of satisfaction** in our patients.

High clinical efficacy results evidenced by **parameters of vaginal health and studies *in vivo***.

Efficacy demonstrated in studies performed by external organisms in **women's intimate rejuvenation**.

Solutions developed by specialist gynaecologists for **function and aesthetic** regeneration of women's intimate area, providing **comfort and improving quality of life**.



1. Gandhi J, Chen A, Dagur G, Suh Y, Smith N, Cali B, Ali Khan S. Genitourinary syndrome of menopause: an overview of clinical manifestations, pathophysiology, etiology, evaluation, and management. Am J Obstet Gynecol. 2016;215(6):704-711. 2. Gloria Bachmann, MD, Richard J Santen, MD. Clinical manifestations and diagnosis of genitourinary syndrome of menopause (vulvovaginal atrophy). Uptodate: Robert L Barbieri, MD. [Internet]. 2019. [Reviewed on 11 February 2020]. 3. Shugusheva L, Gvasalia B, Tkachilova S, Sarakhova D. Efficacy of injectable hyaluronic acid for the treatment of dyspareunia. Journal of Sexual Medicine 2016 13:5 SUPPL. 2 (S121-). 2. Keen MA. Hyaluronic acid in dermatology. Skinmed 2017;15:441-448. 4. Gaspar A, Palikov K, Luis Shotze Luis G. An open-label uncontrolled single centre study for the evaluation of the performance characteristics (efficacy and safety) of mesofiller® sensitive (20 mg/ml) for female intimate area hydration and atrophy of the vagina. 5. Gaspar A, Tsvetkov D, Luis Shotze Luis G. An open-label uncontrolled single centre study for the evaluation of the performance characteristics (efficacy and safety) of mesofiller® sensitive plus (25 mg/ml) female intimate area augmentation of labia majora. 6. Simmler, C., Pauli, G. F., & Chen, S.-N. (2013). Phytochemistry and biological properties of glabridin, 90, 160-184. doi:10.1016/j.fitote.2013.07.003. 7. YOKOTA, T., NISHIO, H., KUBOTA, Y., & MIZOGUCHI, M. (1998). The Inhibitory Effect of Glabridin from Licorice Extracts on Melanogenesis and Inflammation. Pigment Cell Research, 11(6), 355-361. 8. Bisset et al, 2003. 9. Pelletier-Louis, M.-L. (2017). Peelings chimiques et prise en charge du vieillissement cutané. Annales de Chirurgie Plastique Esthétique, 62(5), 520-531. doi:10.1016/j.anplas.2017.07.001. 10. Unit of Biotechnology. Study *in vitro* on melanin synthesis by spectrophotometry mesoestetic Pharma Group, S.L. 2019. 11. Unit of Biotechnology. Study *in vitro* of the efficacy of Glycyrrhiza Glabra root Extract in the inhibition of melanin. mesoestetic Pharma Group, S.L. 2019. 12. Unit of Biotechnology. Study *in vitro* of the efficacy of niacinamide in the prevention of transepidermal water loss TEWL. mesoestetic Pharma Group, S.L. 2019. 13. Unit of Biotechnology. Study *in vitro* of the efficacy of Plankton extract in the inhibition of melanin. mesoestetic Pharma Group, S.L. 2019. 14. Unit of Biotechnology. Study *in vitro* of the efficacy of Plankton extract in the inhibition of melanin. mesoestetic Pharma Group, S.L. 2019. 15. Study *in vivo* of the efficacy of niacinamide in the prevention of transepidermal water loss TEWL. mesoestetic Pharma Group, S.L. 2019. 16. A multicentre, efficacy and safety clinical study of a new method (dermamelan® intimate) consisting of a 2-stage protocol: one session at the clinic and a treatment for home use



mesoestetic Pharma Group, s.l.
C/ Tecnologia, 25
08840 - Viladecans (Barcelona)
Tel. +34 93 325 20 30
www.mesoestetic.com
made in Spain